




P-Mix Dispersing and Mixer for Laboratory

Evolution Series



Designed for saving cost while maintenance capacity

- ❖ Evolution Price
- ❖ Adjustable high speed inverter
- ❖ Easy operation and maintenance

Option

- Homo/Basket mill type
- Multi purpose for Mixing, dispersing and milling

Model	Motor (HP)	Speed (rpm)	Mixing Cap (liters)	Impeller (mm)	Stroke (mm)	A Width (mm)	B Length (mm)	C Height (mm)	Height + Stroke (mm)	Weight (Kg)
Evo 0.1	0.1	0 - 3000	0.2 - 2	35	150	260	255	660	810	14
Evo 0.5	0.5	0 - 3000	0.5 - 4	50	215	515	330	790	1000	21
Evo 1	1	0 - 3000	1.0 - 10	75	260	630	400	920	1160	39
Evo 1-3 Ø	1	0 - 3000	1.0 - 20	100	370	730	450	1060	1300	40

P.General Products : 234 Moo.5, Prakasa Road., Prakasa, Muang, Samutprakarn 10280 Thailand.

Tel: (662) 728-5246-8 Fax: (662) 728-5249

E-mail: info@p-general.com